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








mRNA 	<ul style="list-style-type: none">• Low protein expression• Immunogenicity favorable
Carrier 	<ul style="list-style-type: none">• Protection against RNase• No tropism necessary• Adjuvant immunogenic capacity
In vivo delivery 	<ul style="list-style-type: none">• Local intramuscular delivery• Low mRNA dosage• Low toxicity risks
Clinical target 	<ul style="list-style-type: none">• Immune protection• Circulating antibodies• Immune memory

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mRNA 	<ul style="list-style-type: none">• Prolonged mRNA stability• High protein expression• Minimization of immunogenicity
Carrier 	<ul style="list-style-type: none">• Protection against RNase• Cell- or tissue-specific tropism• Immune cloaking
In vivo delivery 	<ul style="list-style-type: none">• Solid organ direct delivery• Systemic delivery (liver)• High mRNA dosage• High delivery efficiency• Higher toxicity risks
Clinical target 	<ul style="list-style-type: none">• Enzyme replacement• Solid cancer• Tissue regeneration• Gene editing• Antibody therapy

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
Analytical Procedures for mRNA Vaccine Quality (USP Draft Guidelines 2nd Edition)

Identity	mRNA sequence identity confirmation	HIGH-THROUGHPUT SEQUENCING (HTS) Sanger sequencing RT-PCR
Content	RNA concentration	qPCR dPCR UV
Integrity	mRNA intactness	CAPILLARY ELECTROPHORESIS CAPILLARY GEL ELECTROPHORESIS AGAROSE GEL ELECTROPHORESIS
Purity	5' capping efficiency	RP-LC-MS/MS IP-RP-HPLC
	3' poly(A) tail length	IP-RP-HPLC
	Product related impurities - dsRNA	IMMUNOBLOT ELISA
	Product related impurities - aggregate quantitation	SEC-HPLC
	Product related impurities - percentage of fragment mRNA	RP-HPLC
	Process related impurities-residual DNA template	qPCR
	Process related impurities - quantitation of free/ non-incorporated nucleosides	RP-LC-MS/MS
	Process related impurities - residual T7 RNA polymerase content	ELISA
Potency	Expression of target protein	CELL-BASED ASSAY

Quality	Attribute	Method	
Identity	mRNA sequence identity confirmation	High throughput sequencing (HTS)	
		Sanger sequencing	
		Reverse Transcriptase - PCR (RT-PCR)	
Content	RNA concentration	Quantitative PCR (qPCR)	
		Digital PCR (dPCR)	
		Ultraviolet Spectroscopy (UV)	
Integrity	mRNA intactness	Capillary electrophoresis²	
		Capillary gel electrophoresis (CGE)²	
		Agarose gel electrophoresis	
Purity	5' capping efficiency	Reverse-phase liquid chromatography mass spectrometry (RP-LC-MS/MS)²	
		Ion pair reversed-phase high-performance liquid chromatography (IP-RP-HPLC)	
	3' poly(A) tail length	Ion pair reversed-phase high-performance liquid chromatography (IP-RP-HPLC)	
		Immunoblot	
	Product related impurities - dsRNA	Product related impurities - dsRNA	Enzyme-linked immunosorbent assay (ELISA)
			Size exclusion-high-performance liquid chromatography (SEC-HPLC)²
			Reversed-phase HPLC (RP-HPLC)²
			quantitative PCR (qPCR)
Reverse-phase liquid chromatography mass spectrometry (RP-LC-MS/MS)²			
Process related impurities - residual T7 RNA polymerase content	Process related impurities - residual T7 RNA polymerase content	Enzyme-linked immunosorbent assay (ELISA)	
Potency	Expression of target protein	Cell-based assay	
Safety	Endotoxin	USP <85>	
	Bioburden	USP <61>, <62>, <1115>	
Other	Appearance	USP <790>	
	Residual solvents	USP <467>	
	pH	USP <791>	

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
Analytical Procedures for mRNA Vaccine Quality (USP Draft Guidelines 2nd Edition)

Identity	mRNA sequence identity confirmation	Sanger sequencing
		RT-PCR
	Identity of lipids	RP-HPLC-CAD
Content	RNA concentration/RNA encapsulation efficiency	Fluorescence-based assay
	Lipid content	RP-HPLC-CAD
Integrity	LNP size and polydispersity	DYNAMIC LIGHT SCATTERING
	RNA size and integrity	CAPILLARY GEL ELECTROPHORESIS
Purity	Product related impurities - aggregate quantitation	SEC-HPLC
	Product related impurities - percentage of fragment mRNA	IP-RP-HPLC
Potency	Expression of target protein	CELL-BASED ASSAY

Quality	Attribute	Method
Identity	mRNA sequence identity confirmation	Sanger sequencing
		Reverse Transcriptase - PCR (RT-PCR)
	Identity of lipids	Reversed-phase high-performance liquid chromatography with charged aerosol detector (RP-HPLC-CAD)
Content	RNA concentration/RNA encapsulation efficiency	Fluorescence-based assay
	Lipid content	Reversed-phase high-performance liquid chromatography with charged aerosol detector (RP-HPLC-CAD)
Integrity	LNP size and polydispersity	Dynamic light scattering (DLS)
	RNA size and integrity	Capillary gel electrophoresis (CGE)²
Purity	Product related impurities - aggregate quantitation	Size exclusion-high-performance liquid chromatography (SEC-HPLC)²
	Product related impurities - percentage of fragment mRNA	Ion pair reversed-phase high-performance liquid chromatography (IP-RP-HPLC)²
Potency	Expression of target protein	Cell-based assay
Safety	Endotoxin	USP <85>
	Sterility	USP <71>
Other	Appearance	USP <790>
	Residual solvents	USP <467>
	Osmolality	USP <785>
	Subvisible particles	USP <787>
	Residual solvents	USP <467>
	Extractable volume	USP <1>, <698>
	Container closure integrity	USP <1207>
	pH	USP <791>

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