Agilent EZChrom *Elite*

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Agilent Technologies

Notices

Copyright © Scientific Software, Inc 1997-2003 © SSI 2006.

No part of this manual may be reproduced in any form or by any means (including electronic storage and retrieval or translation into a foreign language) without prior agreement and written consent from Agilent Technologies, Inc. as governed by United States and international copyright laws.

Edition

March, 2006

Document Revision 3.2

Printed in USA

Agilent Technologies, Inc.
6612 Owens Dr.
Pleasanton, CA 94588-3334

Warranty

The material contained in this document is provided "as is," and is subject to being changed, without notice, in future editions. Further, to the maximum extent permitted by applicable law, Agilent disclaims all warranties, either express or implied, with regard to this manual and any information contained herein, including but not limited to the implied warranties of merchantability and fitness for a particular purpose. Agilent shall not be liable for errors or for incidental or consequential damages in connection with the furnishing, use, or performance of this document or of any information contained herein. Should Agilent and the user have a separate written agreement with warranty terms covering the material in this document that conflict with these terms, the warranty terms in the separate agreement shall control.

Technology Licenses

The hardware and/or software described in this document are furnished under a license and may be used or copied only in accordance with such license.

Restricted Rights Legend

If software is for use in the performance of a U.S. Government prime contract or subcontract, Software is delivered and licensed as "Commercial computer software" as defined in DFAR 252.227-7014 (June 1995), or as a "commercial item" as defined in FAR 2.101(a) or as "Restricted computer software" as defined in FAR 52.227-19 (June 1987) or any equivalent agency regulation or contract clause. Use, duplication or disclosure of Software is subject to Agilent Technologies' standard commercial license terms, and non-DOD Departments and Agencies of the U.S. Government will receive no greater than Restricted Rights as defined in FAR 52.227-19(c)(1-2) (June 1987). U.S. Government users will receive no greater than Limited Rights as defined in FAR 52.227-14 (June 1987) or DFAR 252.227-7015 (b)(2) (November 1995), as applicable in any technical data.

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```
Dim nTrace
Dim Peaks

nTrace = 1
Set Peaks = CDS.Peaks(nTrace, cNamedDetPeaks)
MsgBox "There are " & Peaks.Count & "named detected peaks"
```

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```
Trace = 1
nPeakID = 2001
Set Peaks = CDS.Peaks(Trace, cNamedDetPeaks)
If(VarType(Peaks) <> vbEmpty) Then
    Set Peak = Peaks.PeakID(nPeakID)
    If(VarType(Peak) <> vbEmpty) Then
        'Do something with the peak
    End If
End If
```

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```
Trace = 1
nPeakIndex = 3
Set Peaks = CDS.Peaks(Trace, cNamedDetPeaks)
If(VarType(Peaks) <> vbEmpty) Then
    Set Peak = Peaks.PeakIndex(nPeakIndex)
    If(VarType(Peak) <> vbEmpty) Then
        'Do something with the peak
    End If
End If
```

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```
Trace = 1
strName = "Peak 1"
Set Peaks = CDS.Peaks(Trace, cNamedDetPeaks)
If(VarType(Peaks) <> vbEmpty) Then
    Set Peak = Peaks.PeakName(strName)
    If(VarType(Peak) <> vbEmpty) Then
        'Do something with the peak
    End If
End If
```

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```
Dim Area
Dim Height
Set Peaks = CDS.Peaks(nTrace, cNamedDetPeaks)
If(VarType(Peaks) <> vbEmpty) Then
  For nPeak = 1 To Peaks.Count
    Peak = vbEmpty
    Set Peak = Peaks.PeakIndex(nPeak)
    If(VarType(Peak) <> vbEmpty) Then
      Peak.GetResult cPeakArea, Area
      Peak.GetResult cPeakHeight, Height
      Peak.SetCustomResult Area/Height
    End If
  Next
End If
```

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```
Dim Area
Dim Height
Set Peaks = CDS.Peaks(nTrace, cNamedDetPeaks)
If(VarType(Peaks) <> vbEmpty) Then
    For nPeak = 1 To Peaks.Count
        Peak = vbEmpty
        Set Peak = Peaks.PeakIndex(nPeak)
        If(VarType(Peak) <> vbEmpty) Then
            Peak.GetResult cPeakArea, Area
            Peak.GetResult cPeakHeight, Height
            Peak.SetCustomResult Area/Height
        End If
    Next
End If
```

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```
Dim nTrace
Dim Peaks

nTrace = 1
Set Groups = CDS.Groups(nTrace,cNamedPeaks)
MsgBox "There are " & Groups.nCount & "named peaks in this group"
```

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```
Trace = 1
Set Groups = CDS.Groups(Trace,cNamedPeaks)
For Count = 0 To Groups.nCount - 1
    Set Group = Groups.Group(Count)
    'Do something with the group object
Next
```

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```
Trace = 1
Set Groups = CDS.Groups(Trace,cNamedPeaks)
For Count = 0 To Groups.nCount - 1
  Set Group = Groups.Group(Count)
  'Do something with the group object
  If(VarType(Group) <> vbEmpty) Then
    Group.GetResult cGroupHeight, Height
    Group.GetResult cGroupArea, Area
    Group.SetCustomResult Area\Height
  End If
Next
```

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```
Trace = 1
Set Groups = CDS.Groups(Trace,cNamedPeaks)
For Count = 0 To Groups.nCount -1
  Set Group = Groups.Group(Count)
  'Do something with the group object
  If(VarType(Group) <> vbEmpty) Then
    Group.GetResult cGroupHeight, Height
    Group.GetResult cGroupArea, Area
    Group.SetCustomResult Area\Height
  End If
Next
```

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```
Set Instrument = CDS.Instrument()  
MsgBox "Employee " & Instrument.UserName & " opened this instrument."
```

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```
Set Project = CDS.Project()  
MsgBox "The data is stored in the folder: " & Project.DataPath
```

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```
Set Project = CDS.Project()  
MsgBox "The project description" & Project.Description
```

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```
Set Project = CDS.Project()  
MsgBox "The methods are stored in the folder: " & _ Project.MethodPath
```

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```
Set Project = CDS.Project()  
MsgBox "Project name: " & Project.Name
```

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```
Set Project = CDS.Project()  
MsgBox "Project root path: " & Project.RootPath
```

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```
Set Project = CDS.Project()  
MsgBox "Project sequence path: " & Project.SequencePath
```

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```
Dim Traces  
Set Traces = CDS.Traces()  
MsgBox "There are " & Traces.Count & "traces."
```

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```
Set Traces = CDS.Traces()  
For Count = 1 To Traces.Count
```

```
Set Trace = Traces.Trace(Count)
'Do something with the Trace object
Next
```

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```
Set Traces = CDS.Traces()
For Count = 1 To Traces.Count
  Set Trace = Traces.Trace(Count)
  MsgBox "The actual point count is " & _      Trace.ActualPointCount
Next
```

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```
Set Traces = CDS.Traces()
For Count = 1 To Traces.Count
  Set Trace = Traces.Trace(Count)
  MsgBox "The actual run time is " & Trace.ActualRunTime
Next
```

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```
Set Traces = CDS.Traces()
For Count = 1 To Traces.Count
    Set Trace = Traces.Trace(Count)
    MsgBox "The base frequency = " & Trace.BaseFrequency
Next
```

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```
Set Traces = CDS.Traces()
For Count = 1 To Traces.Count
    Set Trace = Traces.Trace(Count)
    MsgBox "The Delay time of this trace = " & Trace.DelayTime
Next
```

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```
Set Traces = CDS.Traces()
For Count = 1 To Traces.Count
    Set Trace = Traces.Trace(Count)
    MsgBox "The Detector name this trace =" & Trace.DetectorName
Next
```

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```
Set Traces = CDS.Traces()  
For Count = 1 To Traces.Count  
    Set Trace = Traces.Trace(Count)  
    MsgBox "Expected point count for trace = " & Trace.ExpectedPointCount  
Next
```

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```
MsgBox "Expected run time: " & Trace.ExpectedRunTime
```

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```
StartTime = 12 * 60 ' Time in seconds
StopTime = 18 * 60 ' Time in seconds
Trace.GetMinMaxPoints StartTime, EndTime, MinAmp, MaxAmp,
    vMinAmpTime, vMaxAmpTime
MsgBox "Highest point: " & MaxAmp & " Occurred at: " & vMaxAmpTime
MsgBox "Lowest point: " & MinAmp & " Occurred at: " & vMinAmpTime
```

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GetTracePoints vXPoints, vYPoints

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```
Dim vTime
Trace.IndexToTime 3, vTime
MsgBox "The 3rd point occurred at: " & vTime & " seconds"
```

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```
MsgBox "Sampling period for the trace: " & Trace.SamplingPeriod
```

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```
MsgBox "Trace start time: " & Trace.StartTime
```

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```
MsgBox "Trace stop time: " & Trace.StopTime
```

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```
If (Trace.TraceType = cAuxTrace) Then  
    MsgBox "This is an auxiliary trace"  
Else If  
    MsgBox "This is not an auxiliary trace"  
End If
```

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```
If(Trace.UniformSampling = True) Then  
  MsgBox "Data was uniformly sampled"  
Else  
  MsgBox "Data was not uniformly sampled"  
End If
```

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```
MsgBox "The X axis multiplier is: " & Trace.XAxisMultiplier
```

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```
MsgBox "The X axis title is: " & Trace.XAxisTitle
```

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```
MsgBox "The X axis units is: " & Trace.XAxisUnits
```

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```
MsgBox "The Y axis multiplier is: " & Trace.YAxisMultiplier
```

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```
MsgBox "The X axis title is: " & Trace.YAxisTitle
```

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```
MsgBox "The Y axis units is: " & Trace.YAxisUnits
```

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```
boolean GetSetting(nRequest As Long, vValue)
```

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