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1. Ministry of Health of the People's Republic of China. The Pharmacopeia of the People's Republic of China; 2020 Edition - English Version. <https://store.usp.org/product/3124465>.
2. European Directorate for the Quality of Medicines & HealthCare. European Pharmacopeia. <https://pheur.edqm.eu/home>.
3. Pharmaceuticals and Medical Devices Agency. Japanese Pharmacopeia. <https://www.pmda.go.jp/english/rs-sb-std/standards-development/jp/0029.html>.
4. United States Pharmacopeial Convention. United States Pharmacopeia. <https://www.uspnf.com>.
5. Weber, J.; Buske, J.; Mäder, K.; Garidel, P.; Diederichs, T. Oxidation of Polysorbates: An Underestimated Degradation Pathway? *Int. J. Pharm.X* **2023**, *6*, 100202–100232. DOI: 10.1016/j.ijpx.2023.100202.
6. U.S. Department of Health and Human Services, Food and Drug Administration, Center for Drug Evaluation and Research. Guidance for Industry INDs for Phase 2 and Phase 3 Studies, Chemistry, Manufacturing, and Controls Information; rev. 05.2003, last revised February 2018. <https://www.fda.gov/regulatory-information/search-fda-guidance-documents/inds-phase-2-and-phase-3-studies-chemistry-manufacturing-and-controls-information>.
7. Official Journal of the European Communities. Directive 2001/83/EC of the European Parliament and of the Council of 6 November 2001 on the Community Code Relating to Medicinal Products for Human Use; rev. 11.2001. <https://eur-lex.europa.eu/LexUriServ/LexUriServ.do?uri=OJ:L:2001:311:0067:0128:en:PDF>.
8. Committee for Medicinal Products for Human Use, European Medicines Agency. Information for the Package Leaflet Regarding Polysorbates Used as Excipients in Medicinal Products for Human Use; EMA/CHMP/190743/2016, last revised December 2023. https://www.ema.europa.eu/en/documents/scientific-guideline/information-package-leaflet-regarding-polysorbates-used-excipients-medicinal-products-human-use_en.pdf.
9. Host Cell Protein Contaminants in mAb and Protein Therapy Manufacturing. <https://www.usp.org/biologics/host-cell-proteins>.
10. Aryal, B.; Lehtimäki, M.; Rao, V. A. Stress-Mediated Polysorbate 20 Degradation and Its Potential Impact on Therapeutic Proteins. *Pharm. Res.* **2024**, *41*, 1217–1232. DOI: 10.1007/s11095-024-03700-7.
11. 2023 FDA Science Forum: Polysorbate 20 Degradation in Biotherapeutic Formulations and Its Impact on Protein Quality. <https://www.fda.gov/science-research/fda-science-forum/polysorbate-20-degradation-biotherapeutic-formulations-and-its-impact-protein-quality>.
12. Koppolu, V.; Vemulapalli, B.; Thomas, J.; Wang, S.C.; Borman, J. A Universal Method for the Determination of Polysorbate 80 in Monoclonal Antibodies and Novel Protein Therapeutic Formulations. *Anal. Methods* **2018**, *10*, 1296–1304. DOI: 10.1039/c7ay02537h.
13. Wuchner, K.; Yi, L.; Chery, C.; Nikels, F.; Junge, F.; Crotts, G.; Rinaldi, G.; Starkey, J.; Bechtold-Peters, K.; Shuman, M.; *et al.* Industry Perspective on the Use and Characterization of Polysorbates for Biopharmaceutical Products Part 1: Survey Report on Current State and Common Practices for Handling and Control of Polysorbates. *J. Pharma.Sci.* **2022**, *111*, 1280–1291. DOI: 10.1016/j.xphs.2022.02.009.
14. Wuchner, K.; Yi, L.; Chery, C.; Nikels, F.; Junge, F.; Crotts, G.; Rinaldi, G.; Starkey, J.; Bechtold-Peters, K.; Shuman, M.; *et al.* Industry Perspective on the Use and Characterization of Polysorbates for Biopharmaceutical Products Part 2: Survey Report on Control Strategy Preparing for the Future. *J. Pharma.Sci.* **2022**, *111*, 2955–2967. DOI: 10.1016/j.xphs.2022.08.021.
15. Barnett, A. The Cost of Quality: Can We Really Afford to Ignore It?; *NSF International*: York, UK, **2017**.

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© Agilent Technologies, Inc. 2025